

## Shamir Autograph III®

### Disponibilidad de producto y Marcas láser

Material	Indice	Rangos Esfera	Cilindro	Add	Color	Altura Mínima	Marcas Invisibles
CR-39	1.49	- 6.00 a + 6.00	-5.00	0.75 a 4.00	Blanco	11mm	15 S A <sub>3</sub> ADD
Transitions®	1.49	- 6.00 a + 6.00	-5.00	0.75 a 4.00	Gris & Café	11mm	TR S A <sub>3</sub> ADD
PhotoGraph Colors®	1.56	- 6.00 a + 6.00	- 5.00	0.75 a 4.00	Gris, Café, Rosa, Azul , Morado	11mm	TR S A <sub>3</sub> ADD
Poly Plus	1.59	- 8.00 a + 6.00	- 6.00	0.75 a 4.00	Blanco	11mm	PC S A <sub>3</sub> ADD
Transitions®	1.59	- 8.00 a + 6.00	- 6.00	0.75 a 4.00	Gris & Café	11mm	PT S A <sub>3</sub> ADD
PhotoGraph®	1.59	- 8.00 a + 6.00	- 6.00	0.75 a 4.00	Gris & Café	11mm	PT S A <sub>3</sub> ADD
Polarized	1.59	- 8.00 a + 6.00	- 6.00	0.75 a 4.00	Gris & Café	11mm	PPZ S A <sub>3</sub> ADD
Drivewear®	1.59	- 8.00 a + 6.00	- 6.00	0.75 a 4.00	Drivewear®	11mm	PD S A <sub>3</sub> ADD
Ultra Thin	1.74	- 10.00 a + 6.00	- 7.00	0.75 a 4.00	Blanco	11mm	1.74 S A <sub>3</sub> ADD

\* La curva base real sera definida por el Software Freeform



**SHAMIR**

ReCreating Perfect Vision®

